

THE 'WHY' AND 'HOW' OF A DEDICATED ELECTRON DIFFRACTOMETER:

A user-friendly device for nano-crystallography experiments

Dr. Eric Hovestreydt, Dr. Arianna Lanza, Dr. Gunther Steinfeld, Dr. Gustavo Santiso-Quinones

ELDICO Scientific AG, 5234 Villigen, Switzerland

hovestreydt@eldico.ch | www.eldico-scientific.com

Introduction:

For some reason, no dedicated Electron Diffractometer has been available commercially so far. Data quality would greatly benefit from a setup that focuses on the diffraction capability over imaging and allowing for faster and more complete datasets^[1].

ELDICO's objective has been to create a system by scientists for scientists, paying special attention to the performance, versatility and ease-of-use, which resulted in the smart combination of a 5-axis, 360° rotation, nanometer-precise goniometer and an electron beam of radically simplified design, showing that measurements indeed provide good results and structures can be solved ab-initio.




Unique properties:

- Device optimized for diffraction experiments
- Low dose for beam-sensitive materials
- Stable goniometer for precise nano rotation
- Hybrid-pixel detector for fast acquisition of data


Unique benefits:

- As-easy-to-use-as single-crystal X-ray diffraction
- Straightforward to install – no special requirements
- Sample size from 10 to 1000 nm
- R_1 down to 10 %


First structure obtained on our device: proof of concept



Crystal (STEM/BF)



Diffraction pattern



Structure

$\text{Na}_2\text{Ca}_3\text{Al}_2\text{F}_{14}$

- acceleration: 160kV
- 600 frames: 0.2° / 0.1 sec each
- unit cell: $a = b = c = 10.22(1) \text{ \AA}$
- space group: $I2_13$
- structure solution: Superflip^[2]

[1] M. Gemmi, E. Mugnaioli, T. E. Gorelik, U. Kolb, L. Palatinus, P. Boullay, S. Hovmöller, J. P. Abrahams; ACS Cent. Sci., 2019, 5, 1315-1329.

[2] L. Palatinus, G. Chapuis; J. Appl. Cryst., 2007, 40, 786-79.